

Supore

ACP-800

Auto Chart Projector



NEW

optohellas
Worldwide Representation Platform

Auto Chart Projector

ACP-800

Exquisite appearance and excellent workmanship according to the latest national standard, providing 33charts.

- Fast response in 0.3 seconds
- Setting for red & green LED brightness
- Switching power supply for more stable and reliable machine
- The latest international visual signs for more accurate testing result



Auto Chart Projector

ACP-800

Specification

Chart	33 Selections
Chart selection speed	Average 0.3 sec
Mask	1open
	3horizontal lines
	5vertical lines
	15 single letters
Filter	Red/Green
Projection distance	2.0m-7.0m (5m is standard)
Projection magnification	30X±5X (at 5m)
Projection size	Φ252mm x 330(W) x 225(H)mm (5m)
Tilt angle	Ball joint
Lamp	5V5W, 2000 hrs lifetime
Power source	AC100-240V 50/60Hz
Power consumption	50VA
Auto-Off function	3 minutes after

Dimensions

Main body	265(H) x 230(W) x 360(D)
Remote	20(H) x 61(W) x 180(D)

Net weight

Main body	5.0Kg
Remote	180g

Accessories

	1 dust cover
	1 polarized metal screen
	2 spare fuses
	2 batteries (AAA sized)

Optional accessories

	Table stand
	Wall bracket