

**A.R.C.
LASER**

enlighten your surgery.

**SLT + Nd:YAG
= COBRA**

PRELIMINARY

**A.R.C.
LASER**



COBRA

Nd:YAG & SLT
in ONE SLIT LAMP

SMALL FOOT PRINT
ALL FEATURES OF THE
SINGLE UNITS STAY
TOP TECHNICAL
INTEGRATION



LASER...INNOVATION
MADE IN GERMANY

optohellas®
OPHTHALMIC EQUIPMENT

COBRA



Brilliant Optics for the anterior segment.

5x magnification

Handling via touch screen and dial

SLT with the highest repetition rate worldwide

Precise focus of the Nd:YAG

Made in Germany



Choose SLT or YAG operation via Jog-dial

- Real 2 in 1 laser system
- Less space, more freedom

info@arclaser.de



Details including release dates and retail prices for these product will be announced at a later date.


**A.R.C.
LASER**

COBRA

Nd:YAG + SLT-LASER

**IDEAL IN
CLINIC &
OFFICE**


Capsulotomy
Iridotomy
Glaucoma treatment



VARIO SLT

Amount	Order #	Description
1 x	LSXX	COBRA: Q-switched Nd:YAG Laser 1064 nm, Q-LAS slit lamp PCL5, internal eye safety filter, parallel binocular tube + Frequency doubled, Q-switched Laser 532 nm, up to 10Hz of repetition rate, internal eye safety filter, slit lamp PCL5, parallel binocular tube
	SL11040	Single table (86 cm width) ②
		2 x dust protection cover

ACCESSORIES FOR MORE COMFORT

PCS=1 ZU11038 Elbow rest, 24 mm	PCS=1 ZU11039 Elbow rest, 48 mm
PCS=1 WE01243 Iridotomy lens, Abraham	PCS=1 WE01244 Capsulotomy lens, Abraham
PCS=1 WE01356 SLT Gonio Lens, Latina	PCS=5 WE08956 Hygienic shield

STANDARD DEVICE FEATURES

Hand trigger	Fixed stand with sliders	LED illumination	Parallel binocular tube

EYE SAFETY

AS01033 safety goggle 445/810/980 1064 nm	PCS=1
AS01003 safety goggle 445/514/532/ 980/1064 nm	PCS=1

DEVICE OPTION - ON REQUEST (EXTRA CHARGE)

PCS=4 ME06340 Double castors	BG04008 Foot switch
-----------------------------------	------------------------