



## L-ord

BENEFIT FROM THE  
ADVANTAGE TO PERFORM  
RETINA LASER PROCE-  
DURES\*\* WITH ALMOST  
ANY SLIT LAMP\*

\* Valid for slit lamps similar to HS

\*\* Only in combination with CLASSIC 514



- Swing in and swing out retractable device
- Easy to attach at your slitlamp
- Laser fiber included for CLASSIC 514



# L-ord

Integrate the retina laser CLASSIC 514 to any slit lamp (HS-type)



## Premium eye protection

Unique daylight compensated protection filter is integrated in the optical system

- reliable
- no fading

## Simple handling

While swinging out the eye protection moves out of line of sight. As protective and reliable as a regular filter.

## Scalable spot size

The simple change of the effect spot size from 50  $\mu\text{m}$  up to a maximum of 500  $\mu\text{m}$  eases your daily routine.

## Solid mechanics

The design suits a number of various slit lamps.

To benefit from its usability and flexibility, you easily can take off the L-ord and put it onto another

slit lamp or just put it onto its holder next to the slit lamp to improve your diagnostics.

The rigid construction is designed for a long service life and promises a maintenance-free operation.

Almost  
for any  
slit lamp\*



\* HS-type of  
slit lamps