

HD-20T

Features

HD-20T is specifically designed to accommodate the ophthalmic diagnostic equipment which require additional monitor and computer to display the images taken by the equipment, such as the OCT, optical biometry, specular microscope, etc. It consists of two separate table tops, one is to seat the diagnostic equipment and the other one is for the keyboard and mouse.

SPECS

Unit Measurements:

1120mm×580mm×(680-880)mm

Tabletop size for equipment: 550mm×550mm

Tabletop size for keyboard and mouse: 550×400mm

Weight: 72Kg

Max load: 50Kg

Voltage: 100-220Vac±10%(50/60Hz)

Power: 220VA

