

It's for everyone.



K-0391 Black



K-0391-G White



K-0391-V Violet

Inami Universal Trial Frame K-0391 is designed to maximise examinee's comfort by minimising the frame weight. This highly durable instrument features quick and easy adjustment for both PD and corneal alignment, while providing a sure mounting for trial lenses.



K-0391 Series

Inami Universal Light Weight Trial Frame

PD : 48 ~ 80mm (± 0.8 mm)

Cornea measuring plate : 2mm interval

Frame adjustable angles : 35° (upward 30° $\pm 5^\circ$, downward 5° + 5° - 2°)

Frame adjustable length : 35mm

Nose pad adjustable angles : 360°

Nose pad adjustable length : 22mm ± 1.5 mm

Maximum available lenses : 10 pcs

Axis of lens : 0° $\pm 5^\circ$, 180° $\pm 5^\circ$

Available colour : Black (K-0391), Violet (K-0391-V), White (K-0391-G)