



i-Optik

CM-1900C

Chart Monitor




Distributed By


optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ

Features

- Polarized (optional) chart for binocular test
- 23-inch high contrast and bright display
- Smooth and quiet operation
- Wireless Remote control
- Wide range of test types
- Precise distance adjustment
- Easy mounting
- Random test sequence



Rotating Astigmatic Fan Chart
Fast and Accurate test for Cylinder and Axis



CV-7200 computerized vision tester



CM-1900P
Flat Screen
Monitor Chart

Support data connection
with CV-7200 computerized
vision tester



Specifications

Display	23" LCD (color) High resolution Monitor
Brightness	250cd/m ²
Chart versions	E,C,ABC, 123, Child, special version charts
Chart acuity	6/3-6/150 20/10-20/500 0.04-2.0
Distance	2.0 ~ 7.0m
Power supply	Input: AC100-240V,50/60Hz
Power consumption	≤33W
Dimension	530mm(W)×390mm(H)×120mm(D)
Weight	5.5Kg (Net. Weight)
Remote	IR
Batteries for remote control	3A