

Electrical Table ET01



TABLE ET01



Technical information

Height (without table top)	Min 87 cm, max 92 cm, stroke 25 cm.
Width	82 cm
Depth	38 cm
Weight	22 kg
Lifting capacity	Max 120 kg
Height adjustment mechanism	Electrical
Height adjustment control	1. Up/down switches, sited on the column. 2. Socket for hand/foot control set.
Power supply voltage/frequency	220-240 V~ 50/60 Hz
2 sockets (~230V, max 200W) for power supply to instruments.	

Optional

- V-shaped table top for two instruments.
- Sliding table top for two instruments.
- Chin rest support plate.
- Foot control set for up/down movement of the chair.
- Hand control set for up/down movement of the chair.