

Motorized instrument table MT01, MT02



MT02

MT01



Official Distributor

optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ

Med
instrus

TABLES



MT02

Motorized Instrument table MT02

Technical information

Height	Min 65 cm Max 105 cm Elevating – 40 cm
Elevating speed	25 mm/s
Lifting capacity	100 kg
Weight	38 kg
Height adjustment mechanism	Electrical
Height adjustment control	UP/DOWN switch
Power supply	~115/230V 50/60Hz



MT01

Motorized Instrument table MT01

Technical information

Height	Min 63 cm Max 103 cm Elevating – 40 cm
Elevating speed	30 mm/s
Lifting capacity	50 kg
Weight	18 kg
Height adjustment mechanism	Electrical
Height adjustment control	UP/DOWN switch
Power supply	~115/230V 50/60Hz

