

REMIDIO MICROSCOPE RECORDING DEVICE



MRD or **Microscope Recording Device** enables recording using high resolution optics and artificial intelligence on an iPhone XR.



patient
management
system

wireless
data transfer
and live
broadcast

product specifications

model	mrd 10
maximum video resolution	4K at 60fps
digital zoom	1.1x – 2.5x
focal adjustment	automatic/manual
white balance and exposure control	automatic/manual
smartphone and camera	iPhone XR (removable)
camera resolution	12 MP or more
video output format	.MOV
video flip	horizontal
device weight	1.1kg
dimensions (w x l x h)	128 x 237 x 99 mm
wireless connectivity	Wifi for file transfer

