


optohellas
Worldwide Representation Network



Neocam iX[®]



Neocam iX[®] is designed in Europe by **Optohellas**.


The all new NEOCAM iX, is manufactured by two high innovative ophthalmic manufacturers, SUOER and MedevPlus.

Neocam iX is the 1st eye care system for paediatric & neonatal ophthalmology.

It is a combination of SUOER's SW8000, Worldwide's Best ROP/FA camera with built-in Fluroangiography (FA), and MedevPlus' iXanner 7vn, World's 1st Portable OCT.

By merging these two cutting-edge technologies, Neocam iX offers a complete set of diagnostic capabilities, such as RoP (Color & FA), OCT (line, raster, volume scans), and soon its Artificial Intelligence Software will become the permanent assistant of the users, which will give them a 2nd opinion in diagnosis.


Macula HD




Optic Nerve Cube



FA Fluroangiography




Retina Imaging



Neocam iX[®]





Specification

Technical Data

- Widefield ROP/FA images (140° to 144°)
- OCT scanning images with 7X7mm wide
- 4.3µm micron optical axial resolution (OCT)
- 8mp resolution analysis (ROP)
- 8.1mm world's smallest objective lense (ROP)
- Integrated Fluoroangiography (ROP)
- Artificial Intelligence ready (OCT), (ROP c.s)
- Macula and Optic Nerve OCT exams
- Portability and Connectivity (Dicom)
- Anterior Segment Module (coming soon)
- Lightweight Probes 270gr (OCT), 470gr (ROP)
- 3 Years Warranty (2+1)