



TRIAL LENS CASE OV-266

Set trial lenses, 266 pieces. The Set consists mainly of positive and negative spherical and cylindrical lenses, while there are also special lenses, with which we can measure the pupil distance, the false blindness, and the sensitivity to colors.

TRIAL LENS CASE OV-232

Set trial lenses, 232 pieces is an essential ophthalmic device for every hospital, ophthalmologist and optician.

*** There is available also TLC 48 pieces in Metal/Plastic.



TRIAL LENS FRAME OV-1101

Trial lens frame OV-1101. Stainless steel frame. PD adjustment range of both eyes 48~80mm.

TRIAL LENS FRAME OV-1102

Trial lens frame OV-1102. Titanium frame, particularly resistant and very light (weight 47gr). With a resistant lens retention system and very easy to use. The high-quality metal and plastic construction guarantees you many years of use.

