



#### ■ Specifications

**PD** : 48 ~ 80mm ( $\pm 0.8$ mm)

**Cornea measuring plate** : 2mm interval

**Frame adjustable angles** : 35°(upward 30° $\pm 5^\circ$ , downward 5°+5°-2°)

**Frame adjustable length** : 35mm

**Nose pad adjustable angles** : 360°

**Nose pad adjustable length** : 22mm  $\pm 1.5$ mm

**Maximum available lenses** : 10 pcs

**Axis of lens** : 0° $\pm 5^\circ$ , 180° $\pm 5^\circ$

**Available colour** : Black (K-0391), Violet (K-0391-V), White (K-0391-G)

