

Ophthalmic Unit

EASY PLUS



Distributed By

optohellas
Worldwide Representation Platform

Ophthalmic Unit

EASY PLUS

EASY PLUS has been carefully designed as a modern, space saving Ophthalmic Unit for our the most important Business Partners across the world.

Technical Specifications

- Easy installation at the Doctor's Office with a choice of right - or left - handed setup.
- 2 Instruments sliding table top with Height Adjustment.
- Electrically elevated COMFORT Chair with Swivel $\sim 85^\circ$ and Reclining $\sim 33^\circ$ Systems.
- Auto Ref-Keratometer, Slit Lamp DC Power Supply and adjusting knobs for a Slit Lamp Illumination built-in the table top.
- Aluminum post suitable for installation of both Chart Projector & Phoropter Arm.
- Efficient, Foldable & Revolving LED Lamp.
- Useful Trial Lens Tray with pivotal arm.

Linear PA-3000 or Counterbalanced PA-1000 Phoropter Arm (Options)

