

Ophthalmic Unit

ECO PLUS



Distributed By

optohellas
Worldwide Representation Platform

Ophthalmic Unit

ECO PLUS

Combines modern design with a full range of extras available, and allows for right- or left- handed setup.

Dimensions

	COMFORT CHAIR	PREMIUM CHAIR
Width	1 051 – 1 357mm	1 037 – 1 357mm
Depth	1 414 – 1 522mm	1 373 – 1 481mm
(reclined chair, extended table top)	(1 885mm)	(2 246mm)
Height		1 817mm
Chair working range	493 – 693mm	507 – 707mm
Table top dimensions		880 x 450mm
Table top height (fixed)		878mm
Table top working range (height adjustable)		878 – 978mm





Technical specification

Supply voltage	100 – 240V AC, 50/60Hz
Power	280VA
LED lamp power	9W
Power supply cord	3m with UniSchuko plug
Max. table top load	50kg
Max. chair load	150kg
Max. power of devices supplied from F-type sockets	450VA
Power supply voltage at the table top sockets (SL/AUX1)	2.5 – 24V DC
Supply current of the slit lamp fixator	0-100mA (LED) or 0-250mA (bulb)
Packaging	1 pallet, 120x80x160cm, ca. 225-270kg (gross)

Available colours:



ECO PLUS WITH COMFORT CHAIR



ECO PLUS WITH PREMIUM CHAIR

