

O-Vision[®]
Ophthalmic equipments

OV-101 HFC
PORTABLE FUNDUS CAMERA



MicroClear Imaging Technology

What you see is what you get
Your best partner for fundus
disease diagnosis


Distributed By

optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ



Low weight 450g, movable and portable.

Rechargeable battery provides up to more than 4 hours of continuous operation.



Images real-time display, two times magnification

model	HFC
FOV	45 ΜΟΙΡΕΣ
MINIMIUM PUPIL	2,5MM
REFRACTIVE COMPENSATION	-20D ΈΩΣ +20D
LIGHT SOURCE	WHITE LED/IR
IMAGE RESOLUTION	1920X1080
FOCUS MODE	MANUAL
SCREEN	3,5" ΕΓΧΡΩΜΗ
POWER	≤ 6VA
STORAGE	8GB, MICRO SD CARD
POWER SUPPLY	3,7V LITHIUM BATTERY
INTERFACE	MINI USB/ WIFI
N. WEIGHT	450GR
G. WEIGTH	2,5 KG
SIZE	160MMX90MMX190MM

