



Distributed By



- High Definition and Stable Image.
- 10 Times Magnification and Real-time Display of High Definition Imaging.
- Brand New Optical System and Large Field of View.
- Large Data Service and Various Image Acquisition Mode.
- Three Image Storage Mode:
Images Real-time Display, 80,000 Images Can Be Saved in Micro SD Memory Card.
- USB Connection to Computer.
- Real-time WIFI Connection to Computer, Mobile Terminal and Cloud Transmission.

Technical Parameters

Size: 350px*300px*550px

Net Weight: 420g (Typical)

Resolution Ratio of Image: 1920×1080

Flash Mode: White LED

LCD Screen: 3.5 inch TFT Colored Screen

Picture Format: JPEG(Photograph)

Interface: Micro USB

Image Transmission: Micro USB , Wifi Transmission

Storage: Micro SD card. Maximum to 32G

Power: Li-on rechargeable Battery 3.6V/3400mAh

Forehead Adjustment Range: ≥ 10mm

Illumination Range: ≥ ±30°

Objective lens: 0.42

Numerical Aperture (N) : 0.05

Slit Width: 0mm ~ 12mm

Aperture Diameter: ∅0.2mm、∅1mm、∅5mm、∅12mm

Filter: Red Free, Cobalt blue, Heat Absorption, No filter