

TONOMETER TAT 80-R



R
RODENSTOCK

Distributed By

optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ

TONOMETER TAT 80-R

Type

Applanation Tonometer

Principal

Goldman Tonometry

Specification

Measuring Range	from 0 to 80mmHg in 2mm increments
Accuracy	+0.5mmHg
Diameter of the pneumatic face	3.06mm
Area of Applanation	7.354 mm ²
Probe carrier line separation	43°

Standard Accessories

Calibration Bar

Prism

Mount Base

Tonometer Base Plate

Compatibility Haag Streit Model Slit Lamp

Guarantee One year (against any manufacturing defect)