




Smart, Slim & Easy to use

LCD Chart LCD-1000P

Official Distributor


optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ

LCD Chart LCD-1000P


Smart, Slim & Easy to use

FEAUTURE

- 24 inch high resolution Display Panel with circular polarization.
- Various types of standard charts, polarized charts for Binocular test, ETDRS chart, Contrast test, 3D test, Color Vision test, Maddox test(LED), etc.
- Synchronize with Digital Refractor DR-900.
- Random function , Program mode, Masks, etc.
- 4 kinds of channels for the remote control are available.



Synchronize with Digital Refractor DR-900



LCD-1000P

DR-900

SPECIFICATION

Monitor	Type	24 Inch TFT LED Backlight
	Resolution	1920(H) x 1080(V) Pixels
Functions	Chart	1. Landolt 2. Tumbling E 3. Numbers 4. Letters 5. Child 6. Special
	Masking	Single / Horizontal Line / Vertical Line
	Filter	Red and Green / Circular Polarized Filter
	Special function	Mirror mode / Contrast / White and Black reverse / Random / Crowd
	Working Distance	1.0m~8.0m (Step 0.25m) *
	Control device	Wireless IR remote controller / Digital refractor DR-900
	Screen saver function	Screen off / Chart show (OFF, 5sec, 10sec, 30sec)
	Screen saver waiting time	5min / 10min / 30min
Specifications	Body	573(W) x 363(V) x 46(D) mm / 4.3 kg
	Remote control	64(W) x 195(V) x 20(D) mm / 140 g
	Power supply (Main unit)	DC 12V 3.4A
	Power	100~240V
	Power consumption	70VA
Standard Accessory	Mount	1 pcs
	Glasses	Pola / Red and Green
	Remote control	1 pc
	Battery	2 pcs (AA 1.5V)
	Adaptor	1 pc
	Power Cable	1 pc

Standard Accessory



Remote control



Pola glass



R & G glass

*Some charts will appear incomplete due to difference between screen resolution and image size.

Design and specifications are subject to change without notice.

Manufacturer

Rexxam

Rexxam Co.,Ltd.

Kagawa factory

958,Ikeuchi,Konan-cho,
Takamatsu,Kagawa 761-1494,Japan

Contact

Eye-care Instruments Sales Dept. Tokyo Office

2-4-2 Kandatsukasa-machi, Chiyoda-ku, Tokyo, 101-0048, Japan

TEL:81-3-6262-9471 FAX:81-3-6262-9472

E-mail:eye@rexam.co.jp URL:http://www.rexxam.co.jp URL:http://www.shin-nippon.co.jp

Distributed by

