

TONOMETER AP 151

SHIN-NIPPON
by **Rexxam**



Distributed By

optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ

TONOMETER AP 151

Type

Applanation Tonometer

Principal

Goldman Tonometry

Specification

Measuring Range	from 0 to 80mmHg in 2mm increments
Accuracy	+0.5mmHg
Diameter of the pneumatic face	3.06mm
Area of Applanation	7.354 mm ²
Probe carrier line separation	43°

Standard Accessories

Calibration Bar	
Prism	
Mount Base	
Tonometer Base Plate	
Compatibility	Haag Streit Model Slit Lamp
Guarantee	One year (against any manufacturing defect)