



Digital Camera System

5Mega pixel CMOS camera with
User-friendly filing software

DCS-700


Official Distributor

optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ

DCS-700

FEATURES

- 5Mega pixel CMOS camera
- User-friendly filing software
- Multi images screen
- Full size screen



User-friendly filing software



Both photo and video can be captured and filed in patient database.

Faster data transfer



By USB 3.0 interface, the observation of live image without stress is available.

Specifications

Sensor Type	CMOS color
Sensor Size	1/2"
Valid Pixel Count	2,560 x 1,920 pixel (5Mega pixel)
Power Rating	USB feeding
Power Consumption	1.7 VA
Interface connector	USB 3.0 micro-B
Input/Output connector	8pins connector
Size	135mm(H) x 76mm(W) x 53mm(D) (excluding part of screw)
Weight	Apporx. 540g

Design and specifications are subject to change without notice.

Manufacturer _____

Distributed by _____

Rexxam

Rexxam Co., Ltd.

958 Ikeuchi, Konan, Takamatsu-shi,
Kagawa 761-1494 Japan

Contact _____

Eye-care Instruments Sales Dept. Tokyo Office

2-4-2 Kandatsukasa-machi, Chiyoda-ku, Tokyo, 101-0048, Japan

TEL:81-3-6262-9471 FAX:81-3-6262-9472

E-mail: eye@rexexam.co.jp

URL: http://www.rexxam.co.jp