

TONOMETER SK-T RQ

Digital Applanation Tonometer



sun
Kingdom **SK MED**
For Your Visible Future

Distributed By

optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ

TONOMETER SK-T RQ

Digital Applanation Tonometer

Type

Applanation Tonometer

Principal

Goldman Tonometry

Specification

Measuring Range	from 0 to 80mmHg in 2mm increments
-----------------	------------------------------------

Accuracy	+0.5mmHg
----------	----------

Diameter of the pneumatic face	3.06mm
--------------------------------	--------

Area of Applanation	7.354 mm ²
---------------------	-----------------------

Probe carrier line separation	43°
-------------------------------	-----

Standard Accessories

Calibration Bar	
-----------------	--

Prism	
-------	--

Mount Base	
------------	--

Tonometer Base Plate	
----------------------	--

Compatibility	Haag Streit Model Slit Lamp
---------------	-----------------------------

Guarantee	One year (against any manufacturing defect)
-----------	---