

SUOER

Delta

Ultrasound Scanner




optohellas[®]

Europe & Middle East Representation



B Scan



A Scan



Full Scale UBM

Delta Ultrasound Scanner

B Scan

Frequency	10MHz, Magnetic driven, noiseless
Scanning Mode	Sector Scanning
Magnify	Multi continuous magnification, Real-Time magnification
Resolution	Lateral $\leq$ 0.3mm Verical $\leq$ 0.2mm
Geometry position precision	Lateral $\leq$ 5% Vertical $\leq$ 3%
Depth	60mm
Enhance the part of vitreous body and retina	
Gain of Probe	30dB-105dB
Scanning Angle	53°
Gray Scale	256
False Color	Multi colors
Measurement Type	multigroup distances, perimeters and areas
Image	multiple curves
Postprocessing	processing, Pseudo-color processing curve
Movies	100 images movie, AVI JPG format image output

A Scan

Frequency	10MHz, with LED
Depth	40mm
Precision	$\pm$ 0.05mm
Measurement	Anterior chamber depth, lens thickness, vitreous body lenght, total lenght and average
Eye Mode	Phakic/Aphakic/Dense/Various IOL
IOL Formula	SRK-II, SRK-T,HOFFER-Q, HOLLADAY, BINKHORST-II, HAIGIS
Stat. Calculation	Average and standard deviation
Store	10 Scanning results for each eye


Others

Display Mode	B, B+B, B+A, A
Hint	Preset keyword
Case Search	Multi-keywords
Working Platform	Windows Sysrem
Useridefined report template	

Full Scale UBM

Frequency	50MHz
Scanning Mode	Wide Range Sector, Scanning Mode, Undistorted, Sulcus-to-sulcus
Scanning Range	16mm*9mm;10*6.5mm
Vertical Precision	$\leq$ 40 μ m
Lateral Precision	$\leq$ 40 μ m
Scanning Lines	1024Lines, 15 μ m between each lines
Geometry	Vertical $\leq$ 3%
Location	Lateral $\leq$ 3%
Precision	Non data interpolation, None distortion Imaging
Display Mode	UBM, UBM+A
System	it have a special
Performance	independent 50 μ m ultrasonic amplification system make the anterior segment image clearer
Working Platform	Windows Sysrem

SUOER



SUOER OPHTHALMOLOGY

M. Alexandrou 50, 60133 Katerini
Branch: Evelpidon 61-63
11362 Athens, Greece

Tel: +30 23510 79750 / +30 210 5750572
e-mail: info@opto.gr
[http: www.opto.gr](http://www.opto.gr)

Distributed By

optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ