

SUOER

Ophthalmic Wide Field Imaging System SW-8000P



optohellas®

Europe & Middle East Representation


SUOEP

Ophthalmic Wide Field Imaging System SW-8000P

Technical data

SW-8000 ROP is a Wide Field camera, very light with easy handling.
Is integrating Light Source, Refractive compensation
and High resolution Imaging system


Probe Diameter	8.1mm (smallest probe in the market)
Probe Weight	527.3gr (the lightest in the market)
Probe connectivity	USB3.0 Plug & Play
Wide Field of View	135° to 144°
Light Source	Color Light Source; Fluorescence Angiography
Center field of view	≥ 30 lp/mm
Middle field of view	$(\pm 22.5^\circ) \geq 20$ lp/mm
Edge field of view	$(\pm 45^\circ) \geq 15$ lp/mm
Illumination Source	white LED



optohellas[®]

Europe & Middle East Representation

SUOER



SUOER OPHTHALMOLOGY

M. Alexandrou 50, 60133 Katerini
Branch: Evelpidon 61-63
11362 Athens, Greece

Tel: +30 23510 79750 / +30 210 5750572
e-mail: info@opto.gr
[http: www.opto.gr](http://www.opto.gr)

Distributed By

optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ