

Supore

ACP-1000

Tabletop Refraction System



NEW

optohellas
Worldwide Representation Platform

Tabletop Refraction System

ACP-1000

Most space-saving refraction unit for convenient use compact appearance, sophisticated design. Filter function

- We believe every sqm of your room is valuable
- Space-saving is important
- ACP-1000 focus on improving professional eye testing feeling and value



Tabletop Refraction System

ACP-1000

Specification

Functions & corresponding parameter Chart	33 Selections
Fixing distance	0.2m
Measure distance	5m
Shade	Separate, horizontal, vertical
Filter	Red/Green
Chart control	By Phablet or view tester
Dimensions for chart window	210(W) x 80(H) mm
Lamp-house	LCD or LED
Power source	Main body: AC100-240V, 50/60Hz,
Power	148VA
Sort	Temperature: 10-40 Centigrade
Opposite humidity	30-75%
Atmospheric pressure	860hpa- 1060hpa
Main body dimensions	450 (W), 555 (H), 470 (D) mm, 20kg
Usage life	10 years

