

SUOER

Non-contact Tonometer SW-5000



SUOER



Non-contact Tonometer SW-5000

Technical data

Measuring Range	1mmHg~60mmHg
Measuring Scale	30mmHg, 60mmHg
Measuring Accuracy	1mmHg;
Measuring Distance	11mm
Focus Method	focus points + focus notification
Focus Mode	three-dimensional auto-focus/manual focus/touch screen focus
Interior Light Fixation	Green LED
Stroke of Moving Track	Left-Right: 80mm Forward-Backward: 40mm Up-Down: 20mm
Display	large colored LCD screen
Output	high speed thermal printer


Unique Features

- ① Integrated ORA (Ocular Response Analyzer)
- ② Unique collection of waveform confidence interval data by weight average of three readings, and indicate low confidence interval results
- ③ Manually focus by touching screen
- ④ Non-contact measurements to avoid cross infection
- ⑤ Integrated 24 hours IOP trend analysis system

optohellas®

Europe & Middle East Representation

SUOER



SUOER OPHTHALMOLOGY

M. Alexandrou 50, 60133 Katerini
Branch: Evelpidon 61-63
11362 Athens, Greece

Tel: +30 23510 79750 / +30 210 5750572
e-mail: info@opto.gr
[http: www.opto.gr](http://www.opto.gr)

Distributed By

optohellas[®]
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ