

assistPro

assistPro F Item No.: 2400 0005
assistPro H Item No.: 2400 0009

The perfect assistant's chair for your OR.

Chassis

- High-grade 5-arm aluminum pedestal
- Excellent handling characteristics thanks to smooth-running double castors
- Optionally available with braking double castors
- Supports loads up to 130 kilograms



Seating system


Comfort Seat

- Ergonomically shaped seat and backrest with lumbar support
- Support for correct, back-friendly sitting posture and reduction of tension caused by prolonged sitting
- Durable, easy-care seat cover made of medical-grade artificial leather to meet the highest OR requirements
- Optionally available with electrically conductive upholstery



Flexibility

Variants



Technical specifications

Lifting column

- Height adjustment assistPro F 530–730 mm
- assistPro H 520–720 mm

Chassis

- Five legged chassis ø 585 mm
- Castors ø 65 mm

Comfort seat

- D x W x H (incl. back part) 385 x 425 x 540 mm
- Seat angle adjustment +5°/-5°
- Back angle adjustment +5°/-25°

Weight

- assistPro F 12 kg
- assistPro H 11 kg

Max. load

130 kg

Subject to technical modifications.

- Sterile, comfortable height adjustment via foottriggered gas pressure spring (assistPro F)
- Classic seat height adjustment via manual trigger underneath the seat (assistPro H)
- Flexible height adjustment from 53 to 73 cm (assistPro F) and 52 to 72 cm (assistPro H)
- Optionally available with arm supports (assistPro F)

Standard color

black
Item No.: 3037 0235



Comfort

"PeriLine" Page 105

Components



assistTrend
Item No.: 2500 0002



assistPro F
Item No.: 2400 0005



assistPro H
Item No.: 2400 0009



assistTrend
Armrest pair
incl. traverse

Item No.: 3300 0000



assistPro F
Armrest pair
incl. traverse

Item No.: 3351 0038



easyPro F
Item No.: 2400 0010



easyPro R
Item No.: 2400 0013



Double rimmed castor,
65 mm, with braking function

Item No.: 3351 0040

"PeriLine"



Product Specification

- Without Phthalates
- Biocompatibility:
cytotoxicity (EN ISO 10 993-5) and
epicutaneous test (EN ISO 10 993-10)
- Urine resistance
- Blood resistance
- Resistance to perspiration
- Resistance to spitting
- Light fastness
- Color fastness to rubbing
- Flame retardant

<i>surgiForce</i>	Comfort seat	2600 0004
<i>surgiForce light</i>	Comfort seat	2600 0006
<i>surgiTrend</i>	Comfort seat	2200 0012
<i>surgiLine</i>	Comfort seat	2200 0013
<i>assistPro F</i>		2400 0005
<i>assistPro H</i>		2400 0009
<i>easyPro F</i>		2400 0010
<i>easyPro R</i>		2400 0013
<i>easyLine F</i>		2400 0020
<i>easyLine R</i>		2400 0022

+ Item No.: 3037 0...

optohellas
Worldwide Representation Platform

				
kastanie 015817 Item No.:240	light brown 016037 Item No.:239	creme 015812 Item No.:225	weiß 010881 Item No.:224	stein 015813 Item No.:236
				
saphir 015814 Item No.:231	delft 010766 Item No.:241	antistatisch schwarz Item No.:223	schwarz 011499 Item No.:235	anthracit 016046 Item No.:234
				
azur 015815 Item No.:230	aubergine 016041 Item No.:229	lavendel 016042 Item No.:228	rot 010878 Item No.:238	orange 016040 Item No.:237
				
kiwi 015816 Item No.:233	lime green 016043 Item No.:232	lemon 016039 Item No.:226	safran 010879 Item No.:227	

The colors displayed on your device may differ slightly from the original colors. If you wish to receive original color charts, please contact us at info@ufsk-osys.com.

Upgrade your OR to perfection.



eyeForce



500 XLE



500 XLE comfort



500 XLS



500 ECO



400 XL



motionLine



motionLine-E



surgiForce



surgiForce light



surgiTrend



surgiLine



stayClean



stayClean roll



stayClean kick bucket



supportLine



assistTrend



assistPro F



assistPro H



easyPro F



easyPro R



easyLine F



easyLine R



assistLine