

surgiTrend



Item No. : 2600 0012 with Comfort seat (Illustration)
Item No. : 2600 0010 with Dynamic seat

For back protection and maximum work flexibility in confined spaces. Sophisticated, space-saving drive mechanism with reduced size for ample foot space combines tilt stability with easy and smooth maneuvering. Outstanding sitting comfort and a wide range of adjusting functions allow every surgeon to find his or her personal best and healthiest working position for any work scenario. Available in two seating system types: Comfort or Dynamic.



Seating systems



Comfort

Seat angle adjustment



Back angle adjustment



Back height adjustment



Upholstery



Comfort

"Periline"
Page 97



Dynamic

Dynamic-function



Lumbar support
angle adjustment



Lumbar support
height adjustment



Dynamic



"Stamskin top"
Page 96

SURGICAL CHAIRS & TABLES SURGEON'S CHAIRS / surgiTrend

Official Distributor

optohellas®
ΟΦΘΑΛΜΟΛΟΓΙΚΑ ΜΗΧΑΝΗΜΑΤΑ

Armrests

- Stable, ergonomically shaped armrests with ample support surface (350 mm x 120 mm)
- Height-adjustable and equipped with a ball joint for 360-degree rotation
- Designed for personalized adjustment to the surgeon
- Optionally available with linear guide and one-hand release



Chassis

- 4 high-quality, smooth-running double castors
- 5th castors for added resistance to tipping
- Optimized foot space
- Central braking system effective on 4 castors
- Integrated foot keypad panel for seat-height adjustment
- Power supply-independent battery box with emergency stop button



Technical specifications

Lifting column

- Drive 3.000 NIPXmm
- Lifting height 200 mm

Chassis

- W x L incl. braking system / foot keypad 600 x 650 mm
- Castors ø 100 mm
- 5th castors ø 65 mm
- Central brake acting on all four castors
- Emergency stop button

Battery unit

- Battery box with 2 batteries (lead-gel) 2 x 12 V, 7,2 Ah (pollution free, maintenance-free, cadmium-free)
- Universal battery charger: Input 100 – 240 V Output 24 V

Seating system: Comfort (Dynamic)

- Min. sitting height 550 mm (600 mm)
- Max. sitting height 750 mm (800 mm)
- Seat angle adjustments +5° / -5° (+12°)
- Back angle adjustments +5° / -25° (-16°)

Weight

45 kg

Max. load

150 kg

Subject to technical modifications.

Technical dimensions

Seating system Comfort



Seating system Dynamic



Accessories



surgiTrend
Dynamic seat Item No. : 2200 0010



surgiTrend
Comfort seat Item No. : 2200 0012



surgiLine
Formed seat Item No. : 2200 0008



surgiLine
Comfort seat Item No. : 2200 0013



Backup battery

Item No.: 3047 0000



surgiTrend
Footswitch plate adapter

Item No.: 3351 1011



surgiLine
Footswitch plate adapter

Item No.: 3351 1014



Arm rest pair
with linear guide and
one-hand release,
factory-fitted

Item No.: 33511 006



Arm rest pair
with one-hand release,
factory-fitted

Item No.: 3351 1004

