

Ophthalmic PocketScope™

Welch Allyn 2.5 V Ophthalmic PocketScopes are lightweight, compact diagnostic instruments for those who want high performance in a slightly smaller package. Welch Allyn PocketScopes use halogen illumination and have many other features found on 3.5 V sets. They are smaller in size but not in utility.




Ophthalmic CompacSet™

This portable ophthalmic instrument is specifically designed to make examining patients easier. The Welch Allyn CompacSet features a retinoscope and an ophthalmoscope built into one compact unit.



#93300 Ophthalmic CompacSet with AA batteries and case



POCKETSCOPE SET

#94830 (94805* export) Ophthalmic PocketScope Set with rechargeable 2.5 V batteries, charging stand, and case



INDIVIDUAL POCKETSCOPE

#12820 PocketScope Ophthalmoscope with AA batteries

#16220 PocketScope Retinoscope with AA batteries

#14220 PocketScope Transilluminator with AA batteries

ACCESSORIES

#03900-U PocketScope Ophthalmoscope lamp

#03400-U PocketScope Transilluminator lamp

#08400-U PocketScope Retinoscope lamp

#72600 Rechargeable 2.5 V PocketScope NiCad battery

#72610 Rechargeable 2.5 V Ophthalmic CompacSet NiCad battery

*Available for export. Please specify country and voltage.
Export ordering example: #94805 – Germany – 220 V